WORKER²

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FACHHANDEL

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T. +49 (0)7021 738 30-0
info@perpedes.de
www.perpedes.com

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For the visually impaired, this document is available in PDF format with this QR code.

CONTENTS

Introduction	26
Production instruction/Type-approved treatment system and its intended purpose	27
Indications/contraindications/WORKER ² variants / Contraindications WORKER ² ELTEN + WORKER ² HAIX® + WORKER ² ATLAS + WORKER ² STEITZ SECURA variants.....	28
Instructions for use and information on risks	29
Cleaning and care/storage/reuse.....	31
Technical specifications and parameters/ Processing instructions WORKER ² + WORKER ² ELTEN + WORKER ² HAIX® + WORKER ² ATLAS + WORKER ² STEITZ SECURA variants	32
Materials and adhesives for alterations to WORKER ² prefabricated foot orthoses	34
Approved adhesive/Assembly and adaptation regulations/ Changes to the orthotic shell	35
Labeling of the custom-made device	36
Changes to the orthotic shell	37
EVA Wedge (conductive).....	38
Supination or pronation wedge.....	40
Disposal/Guarantee	42

INTRODUCTION

Dear clients,

Please read the instructions for use carefully before using the product. Follow all instructions, especially the safety instructions.

Perpedes GmbH
Haertwasen 8-14
73252 Lenningen
Germany

T. +49 (0)7021 738 30-0
info@perpedes.de
www.perpedes.com

PRODUCTION INSTRUCTION

For the incorporation of foot orthoses into safety footwear in accordance with European standards DIN EN ISO 20344 / 20345 and 61340.

These production instructions apply for the WORKER² prefabricated foot orthotic blanks for occupational safety footwear manufactured by Perpedes. The individual work steps are defined and the materials and components to be used for the foot orthoses are specified here. The instructions are required in order to ensure conformity with the EC type examination certificate for occupational safety footwear even after insertion of the foot orthoses. Deviation from these production instructions shall render the EC type examination certificate null and void!

TYPE-APPROVED TREATMENT SYSTEM AND ITS INTENDED PURPOSE

The WORKER² foot orthoses were tested and certified in combination with occupational safety footwear from a variety of manufacturers in accordance with DIN EN ISO 20344 / 20345 and 61340 by recognized test institutes.

The combination of the foot orthoses and ESD occupational safety footwear is therefore type-approved and can be used in static and, depending on the shoe model, even in ESD areas.

The testing applies only for the relevant shoe models from the cooperating shoe manufacturers. The currently approved models of occupational safety footwear can be found at

WORKER²: [www.perpedes.de/downloads/worker²](http://www.perpedes.de/downloads/worker2)

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WORKER² ATLAS: www.perpedes.de/downloads/atlas

WORKER² STEITZ SECURA: www.perpedes.de/downloads/steitz_secura

For orthopedic technicians / orthopedic shoe technicians and other qualified specialists, this results in the advantage of being able to supply in compliance with the requirements of the EU guidelines without any further application for testing procedures (subject to possible variations owing to the legal requirements in the specific country of usage).

INDICATIONS

The WORKER² foot orthoses are suitable for the conservative treatment of foot, knee, hip and lower back problems. Alongside the treatment of foot problems, they can also be used for the prevention of overloading of the entire postural and musculoskeletal system (knee, hip and back area). The treating physician is responsible for the exact diagnosis.

CONTRAINDICATIONS

- Diabetic foot syndrome

WORKER² VARIANTS

FOOT ORTHOSE	ITOM NUMBER		
	completely glued, with transfer print	completely glued, without transfer print	non glued
WORKER ² PRO	WOPRKA_	WOPRKB_	WOPRLA_
WORKER ² CALCA	WOCAKA_	WOCAKB_	WOCALA_
WORKER ² CPX	WOCXKA_	WOCXKB_	WOCXLA_
WORKER ² RIGIDUS	WORIKA_	WORIKB_	WORILA_

WORKER² ELTEN VARIANTS

FOOT ORTHOSE	ITOM NUMBER
	completely glued
WORKER ² PRO ELTEN	WOPRKC_
WORKER ² CALCA ELTEN	WOCAKC_
WORKER ² CPX ELTEN	WOCXKC_

WORKER² HAIX® VARIANTS

FOOT ORTHOSE	ITOM NUMBER
	completely glued
WORKER ² PRO HAIX®	WOPRKD_
WORKER ² CALCA HAIX®	WOCAKD_
WORKER ² CPX HAIX®	WOCXKD_

WORKER² ATLAS VARIANTS

FOOT ORTHOSE	ITOM NUMBER
	completely glued
WORKER ² PRO ATLAS	WOPRKF_
WORKER ² CALCA ATLAS	WOCAKF_
WORKER ² CPX ATLAS	WOCXKF_

WORKER² STEITZ SECURA VARIANTS

FOOT ORTHOSE	ITOM NUMBER
	completely glued
WORKER ² PRO STEITZ SECURA	WOPRKE_
WORKER ² CALCA STEITZ SECURA	WOCAKE_
WORKER ² CPX STEITZ SECURA	WOCXKE_

INSTRUCTIONS FOR USE AND INFORMATION ON RISKS

- The WORKER² prefabricated foot orthoses must be properly customized for the foot and adjusted in the shoe by a qualified specialist.
- Foot orthoses should be worn with (crease-free) socks and should always be used as a pair! Exceptions must be stated by the treating physician.
- The shoes must ensure a sufficient interior volume (shoe length / width, toe space) for the foot as well as for the chosen WORKER² prefabricated foot orthoses. The inlays included in the scope of delivery for the occupational safety footwear must be removed before the WORKER² is inserted.
- The chosen prefabricated foot orthoses need to be cut to match the width and length. The forefoot, mid foot and rear foot must have sufficient flat support from the foot orthoses and should not be on the edge as this can lead to the formation of blisters or calluses on the soles of the feet.
- Clear identification of the left and right foot orthoses must be ensured for the patient / user.
- The foot orthoses should not be exposed to temperatures above 60 °C (e.g. heating, driers, direct sunlight).
- The foot orthoses are not suitable for direct contact with acids and bases. The manufacturer provides no warranty in this case.
- Materials which are used for the foot orthoses are subject to aging and their properties may change. If the coverings, cushioning material or orthotic shells have become damaged, holey or wrinkled as a result of mechanical stress or are broken then the foot orthoses should no longer be worn!

- Foot orthoses are affected by body weight, pressure, torsion and rubbing and therefore have a load-dependent lifespan. Temperatures, sock materials, orthotic materials, footwear and environmental factors can also limit the functional life.
- Checking of the function and condition of the foot orthoses by an orthopedic technician / orthopedic shoe technician or other qualified specialist after a maximum of 6 months is recommended.

CLEANING AND CARE

- Foot orthoses must be removed from the shoes daily in order for moisture to evaporate.
- The surface materials can be cleaned as required with mild washing lotion and a damp sponge or cloth. Soap residues should be removed. Do not use any external heat sources to dry, simply leave to air-dry. Note that the foot orthoses are not machine-washable!
- Mild disinfectant (suitable for human skin) can be sprayed on / used in small amounts. Discoloration of the covering material as a result cannot be excluded and is not grounds for a warranty claim.

STORAGE

- Store in the sealed original packaging between 10 °C and 25 °C without mechanical influences and avoiding direct UV radiation.
- Regular checking of the storage and conditions is recommended.

REUSE

- The product is used for an individually fitted treatment for an individual foot orthoses wearer and should not be reused or subsequently used by other third persons.
- The manufacturer accepts no liability in the event of improper use.

TECHNICAL SPECIFICATIONS / PARAMETERS

Foot orthotic materials may vary depending on the model. However, a foot orthoses is composed of the following main components:

- Thermoplastic deformable, self-supporting orthotic shell
- Padding materials (optional)
- Reinforcement materials (optional)
- Supporting components, e.g. met pad (optional)
- Top cover and base cover

PROCESSING INSTRUCTIONS

The production instructions apply for the use of prefabricated foot orthoses with the following item numbers:

WORKER² VARIANTS

FOOT ORTHOSE	ITOM NUMBER		
	completely glued, with transfer print	completely glued, without transfer print	non glued
WORKER ² PRO	WOPRKA_	WOPRKB_	WOPRLA_
WORKER ² CALCA	WOCAKA_	WOCAKB_	WOCALA_
WORKER ² CPX	WOCXKA_	WOCXKB_	WOCXLA_
WORKER ² RIGIDUS	WORIKA_	WORIKB_	WORILA_

WORKER² ELTEN VARIANTS

FOOT ORTHOSE	ITOM NUMBER
	completely glued
WORKER ² PRO ELTEN	WOPRKC_
WORKER ² CALCA ELTEN	WOCAKC_
WORKER ² CPX ELTEN	WOCXKC_

WORKER² HAIX[®] VARIANTS

FOOT ORTHOSE	ITOM NUMBER
	completely glued
WORKER ² PRO HAIX [®]	WOPRKD_
WORKER ² CALCA HAIX [®]	WOCAKD_
WORKER ² CPX HAIX [®]	WOCXKD_

WORKER² ATLAS VARIANTS

FOOT ORTHOSE	ITOM NUMBER
	completely glued
WORKER ² PRO ATLAS	WOPRKF_
WORKER ² CALCA ATLAS	WOCAKF_
WORKER ² CPX ATLAS	WOCXKF_

WORKER² STEITZ SECURA VARIANTS

FOOT ORTHOSE	ITOM NUMBER
	completely glued
WORKER ² PRO STEITZ SECURA	WOPRKE_
WORKER ² CALCA STEITZ SECURA	WOCAKE_
WORKER ² CPX STEITZ SECURA	WOCXKE_

MATERIALS AND ADHESIVES FOR ALTERATIONS TO WORKER² PREFABRICATED FOOT ORTHOSES

The foot orthoses must only be modified using the adhesives and covering materials specified by the Perpedes GmbH in order to guarantee the height and continuing ESD conductivity.

	6 MM EVA WEDGE (CONDUCTIVE)	10 MM EVA WEDGE (CONDUCTIVE)
Size 1 (35 bis 39)	WOK61	WOK101
Size 2 (40 bis 44)	WOK62	WOK102
Size 3 (45 bis 48)	WOK63	WOK103
Antistat adhesive Ultraflex, 1 kg	79AK1	79AK1
Thinner for Antistat, 1 liter	79IS1	79IS1

Ordering address:

Perpedes GmbH
Härtwasen 8-14, 73252 Lenningen, Germany
T. +49 (0)7021 738 30-0, F. +49 (0)7021 738 30-298
info@perpedes.de, www.perpedes.com

APPROVED ADHESIVE

- Renol AL, Renia GmbH
- Perpedes Antistat adhesive Ultraflex - Item No: 79AK1

ASSEMBLY AND ADAPTATION REGULATIONS

The foot orthoses may only be modified or adapted in accordance with the procedures set out below and as defined.

You can build up and adjust the inlay to meet the requirements of the wearer as long as you take account of the area beneath the toe cap which cannot be altered (up to 10 mm behind the edge of the toe cap - the total thickness of the foot orthoses must not be changed in the area of the steel toe cap).

The largest and flattest possible contact surface with the inner shoe sole / template of the shoe must be ensured through proper grinding and modification of the foot orthotic's width and length in order to fulfill the contact resistance requirements in anti-static and ESD areas. Thanks to the special material combination of the WORKER² foot orthoses, it is possible to individually modify the construction. Functional zones can be built up.

WARNING regarding the HAIX[®] CO system:

The antistatic properties and ESD capability of all EG prototype-tested HAIX[®] professional and safety footwear are secured via an electrically conductible textile strip on the inside of the shoe.

This grey anti-static band must extend at least 5 centimetres into the shoe. In order to continue to be able to guarantee conductivity, the inlay must not cover the strip and the contact

area to the foot needs to be about 2 centimetres. In light of the above information regarding the conductive strip, the use of antistatic or ESD-capable materials (such as EVA / adhesive) is not absolutely necessary.

LABELING OF THE CUSTOM-MADE DEVICE

The orthopedic change must be marked in the shoe below the orthopedic foot orthoses by the specialist retailer:

- a) with the name of the manufacturer
(name of the orthopedic expert)
- b) with the type number

The manufacturer of the custom-made device must apply the procedure for custom-made devices according to Annex XIII of the Regulation (EU) 2017/745 on medical devices (MDR).

CHANGES TO THE ORTHOTIC SHELL

Fully glued foot orthoses:

- by means of a heat gun from the bottom
- approx. 140 °C
- approx. 3 mins
- for surface heating at a distance of 30 cm
- core can be moulded by hand or positively moulded

Heat-based material shrinkage of the EVA bottom cover cannot be excluded. We recommend careful heating with a hot air gun.

Foot orthoses with loose top cover:

- by means of a heat gun from the top side
- approx. 140 °C
- approx. 3 mins
- for surface heating at a distance of 20 cm
- core can be moulded by hand or positively moulded

Heat-based material shrinkage of the EVA bottom cover cannot be excluded. We recommend careful heating with a hot air gun.

EVA WEDGE (CONDUCTIVE)

The following modifications and alterations are possible with the defined EVA materials in accordance with DIN EN ISO 20345 and 61340:

Compensation for leg length discrepancies

For compensating of leg length differences of a max. 10 mm or for bilateral heel elevation of up to 10 mm (Fig. 1). The compensation lift must be carried out proximally to the MTP-joints and may not exceed the original height of the WORKER² prefabricated foot orthoses in the steel toecap area of the occupational safety footwear (Fig. 2)!

Fig. 1

Fig. 2

Processing

1. Cut the conductive EVA wedge (optionally available) to the required size and roughen it on a grinding belt (granulation of 24-40). Thorough de-dusting of the surface to be bonded is then required.
2. Apply an even layer of 79AK1 antistat adhesive or Renol AL to both of the materials to be bonded. Allow the adhesive to dry off for at least 5 minutes.
3. Activate the adhesive film on the build-up material using a heat gun at a temperature of 120 °C for a max. of 20 seconds and on the EVA Worker bottom covering for max. 5 seconds. Then press the build-up material onto the pre-fabricated foot orthoses by hand or with a suitable positive using the vacuum method until it has cooled.

SUPINATION OR PRONATION WEDGE

To modify the WORKER² prefabricated foot orthoses, you can apply and alter the appropriate conductive EVA wedge (optionally available), respecting a maximum height of 4 mm in the heel area, beginning from the lateral or medial heel edge to a maximum of the MTP-joints in a wedge shape to the center line of the foot orthoses at 0 mm (Fig. 3-5). The original height of the WORKER² prefabricated foot orthoses may not be exceeded in the steel toecap area of the occupational safety footwear!

Fig. 3

Fig. 4

Fig. 5

Processing

1. Cut the conductive EVA wedge (optionally available) to the required size and roughen it on a grinding belt (granulation of 24-40). Thorough de-dusting of the surface to be bonded is then required.
2. Apply an even layer of 79AK1 antistat adhesive or Renol AL to both of the materials to be bonded. Allow the adhesive to dry off for at least 5 minutes.
3. Activate the adhesive film on the build-up material using a heat gun at a temperature of 120 °C for a max. of 20 seconds and on the EVA Worker bottom covering for max. 5 seconds. Then press the build-up material onto the pre-fabricated foot orthoses by hand or with a suitable positive using the vacuum method until it has cooled.

DISPOSAL

You may dispose of the product in accordance with the applicable national regulations.

GUARANTEE

The current statutory provisions apply. Disregard of the instructions for use, improperly performed cleaning as well as independent technical changes and additions to the product will void the warranty and the general product liability.

The CE marking confirms compliance with the requirements of Regulation (EU) 2017/745 on medical devices (MDR).

The MDR certificate can be downloaded from our homepage.

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