

# 2-WEITENSYSTEM / 2-WIDTHS SYSTEM

## GRÖSSENBESTIMMUNG OT / OS / AV / NT

## SIZE DETERMINATION OT / FO / AV / NT

Schuhgröße Perpedes Shoe size Perpedes	Fußlänge in cm Foot length in cm	Innensohlen- / Brandsohlenlänge in cm Inner shoe length in cm	Schuhgröße UK UK size
18	10,5	12,1	2K
19	11,2	12,7	3K
20	11,8	13,4	4K
21	12,5	14,0	4,5K
22	13,2	14,7	5,5K
23	13,8	15,3	6K
24	14,5	16,0	7K
25	15,2	16,6	7,5K
26	15,8	17,3	8,5K
27	16,5	18,0	9K
28	17,2	18,7	10K
29	17,8	19,3	11K
30	18,5	20,0	12K
31	19,2	20,7	12,5K
32	19,8	21,4	13,5K
33	20,5	22,0	1
34	21,2	22,7	2
35	21,8	23,4	3
36	22,5	23,7	3,5
37	23,2	24,4	4
38	23,8	25,1	5
39	24,5	25,7	6
40	25,2	26,4	6,5
41	25,8	27,0	7
42	26,5	27,7	8
43	27,2	28,4	9
44	27,8	29,1	9,5
45	28,5	29,8	10,5

Schuhgröße Shoe size	*AV, OS, NT-Fersenmaß *AV, FO, NT heel girth	*AV, OS, NT-Ballenmaß *AV, FO, NT ball girth	*OT-Fersenmaß *OT heel girth	*OT-Ballenmaß *OT ball girth
18	20,2	15,0	22,3	17,5
19	20,8	15,3	22,8	17,8
20	21,3	15,5	23,3	18,0
21	21,8	15,8	23,8	18,2
22	22,4	16,0	24,3	18,5
23	22,9	16,3	24,8	18,7
24	23,5	16,5	25,3	18,9
25	24,0	16,8	25,8	19,2
26	24,7	17,1	26,4	19,5
27	25,5	17,4	27,1	19,9
28	26,2	17,8	27,7	20,3
29	27,0	18,2	28,6	20,7
30	27,7	18,7	29,4	21,2
31	28,5	19,1	30,3	21,6
32	29,2	19,6	31,1	22,1
33	30,0	20,1	32,0	22,6
34	30,6	20,6	32,7	23,1
35	31,1	21,1	33,3	23,6
36	31,7	21,7	34,0	24,2
37	32,2	22,2	34,6	24,7
38	32,8	22,6	35,6	25,1
39	33,5	23,0	36,5	25,5
40	34,1	23,4	37,5	25,9
41	34,8	23,9	38,5	26,4
42	35,4	24,3	39,2	26,8
43	36,0	24,7	39,9	27,2
44	36,7	25,1	40,7	27,6
45	37,3	25,5	41,4	28,0

Für Halbschuhe, Halbstiefel (cm)  
For sneakers, ankle boots (cm)

Schuhgröße Shoe size	*AV, OS, NT-Fersenmaß *AV, FO, NT heel girth	*AV, OS, NT-Ballenmaß *AV, FO, NT ball girth	*OT-Fersenmaß *OT heel girth	*OT-Ballenmaß *OT ball girth
18	19,3	14,5	22,3	17,5
19	19,9	14,8	22,8	17,8
20	20,5	15,0	23,3	18,0
21	21,1	15,2	23,8	18,2
22	21,7	15,5	24,3	18,5
23	22,3	15,7	24,8	18,7
24	22,9	16,0	25,3	18,9
25	23,5	16,2	25,8	19,2
26	24,2	16,6	26,4	19,5
27	24,9	17,0	27,1	19,9
28	25,6	17,4	27,7	20,3
29	26,3	17,8	28,6	20,7
30	27,0	18,2	29,4	21,2
31	27,8	18,5	30,3	21,6
32	28,5	19,9	31,1	22,1
33	29,2	19,3	32,0	22,6
34	29,9	19,9	32,7	23,1
35	30,6	20,4	33,3	23,6
36	31,3	21,0	34,0	24,2
37	32,0	21,5	34,6	24,7
38	32,7	22,1	35,6	25,1
39	33,4	22,6	36,5	25,5
40	34,1	23,2	37,5	25,9
41			38,5	26,4
42			39,2	26,8
43			39,9	27,2
44			40,7	27,6
45			41,3	28,0

Für Sandalen (cm)  
For sandals (cm)



In Therapieschuhen von Perpedes ist eine Spitzenzugabe von 1,2 - 1,5 cm bereits enthalten. Sie messen nur die Fußlänge und lesen die dazugehörige Schuhgröße ab.

Brandsohlenlänge / korrekte Schuhgröße = Fußlänge in cm + Spitzenzugaben (1,2 - 1,5 cm)

Beispiel:  
Gemessene Fußlänge 21,8 cm + Spitzenzugabe = 23,4 cm. Unsere Empfehlung: Größe 35

Perpedes' therapy shoes already include a length allowance of 1.2 - 1.5 cm. You only measure the foot length and read the corresponding shoe size in the table above.

Proper inner shoe length / proper shoe size = foot length in cm + allowance (1.2 - 1.5 cm)

Example:  
Measured foot length of 21.8 cm + allowance = 23.4 cm. Our recommendation: size 35

- \*OT/OT = Orthesenschuhe / AFOs/DAFOs
- \*OS/FO = Orthesenschuhe schmal / Foot orthotics
- \*AV/AV = Antivarusschuhe / Antivarus shoes
- \*NT/NT = Neutralschuhe / Neutral shoes